




Increase department throughput
Drive improved efficiency



Increase in image resolution
Improve diagnostic confidence



Patient comfort in the MR
Faster exams = less table time



Schedule flexibility
Administrative support



Increase clinical referrals
Expand services

Deep Resolve and Value across the Radiology Network

siemens-healthineers.us/deep-resolve

Deep Resolve is the Siemens Healthineers AI-powered image reconstruction technology that takes advantage of convolutional neural networks to accelerate MRI scans, faster than ever before while simultaneously increasing resolution.





Deep Resolve Gain

What is it?

Deep Resolve Gain is a targeted denoising method to increase the MRI signal of images. With this, either shorter scan times or higher resolution can be achieved.

How does it work?

The acquired MRI system generates a patient specific, targeted noise map which reflects spatial noise variations. The MR image and the corresponding targeted noise map are subtracted in an iterative process to remove excess noise which does not add value to the final MRI image.

What problem does it solve?

MRI coil geometry and patient variability can create unnecessary noise in MR images. The noise impairs the signal of the MRI image and can decrease image quality. Deep Resolve Gain directly addresses the removal of the signal not contributing to the image.

What is the key differentiator?

Deep Resolve Gain targets the local noise variation introduced by each MRI patient and the corresponding RF coil via the iterative principle, similar to compressed sensing, but focused on noise removal.



The benefits:



Speed up scan time



Get higher image quality

Deep Resolve Sharp

What is it?

Trained on over 10,000 pairs of low-res and high-res MR data, Deep Resolve Sharp reconstructs a higher resolution image from low resolution data via the AI network.

How does it work?

A low-res image is reconstructed from the low-res data into a sharper, higher resolution image. Data consistency with the acquired raw data is ensured as part of the final image reconstruction. The result is an image with sharper edges and up to a 2x improvement of in-plane resolution.

What problem does it solve?

Typically, in MRI, you can either have a longer scan that acquires a lot of data and provides high resolution or a quicker scan, less data – but lower resolution. Why not have the best of both scenarios with Deep Resolve Sharp.

What is the key differentiator?

The AI network can generate high-res output from low-res input because it has been trained on many pairs of high and low-res images, enabling the software to anticipate where to expect a sharp edge in an image.



The benefits:



Speed up scan time



Get higher image sharpness

especially on the edges

Deep Resolve Boost

What is it?

Deep Resolve Boost delivers super-resolution head to toe.

How does it work?

Leveraging raw data from a reduced, and thus faster scan as input, a deep neural network is applied multiple times in an iterative process to generate the final output with significantly reduced noise. The integration of the raw data along the entire reconstruction process leads to unmatched performance and ensures data integrity.

What problem does it solve?

With conventional reconstruction methods, a highly accelerated image acquisition will lead to strong noise contamination and/or artifacts. The reconstruction with Deep Resolve Boost enables the generation of images with extremely high MRI signal and a superfast image acquisition simultaneously.

What is the key differentiator? Deep Resolve Boost can be combined with Deep Resolve Sharp as well as Simultaneous Multi-Slice (SMS) for amazingly fast accelerations and super-resolution across all anatomies.



The benefits:



Fast acquisition



Super-resolution




Scanner Efficiency with Deep Resolve

Customer Need:

To drive increased throughput while maintaining high image quality

Siemens Healthineers Solution:

Deep Resolve: Acquiring sharper images, faster



MR Exam (1.5T Altea)	Pre Deep Resolve (minute)	Post Deep Resolve (minute)	Exam time savings	% exam time reduction
Brain	12:42	6:27	6:15	49%
T1 Sag	2:17	1:25	:52	38%
DWI Ax	1:14	0:33	:41	55%
T2 Ax	2:10	1:24	:46	35%
T2 IR Ax	4:14	1:30	2:44	65%
T1 Ax	2:47	1:35	1:12	43%
Knee	19:03	8:33	10:30	55%
T2 FS Ax	3:32	1:46	1:46	50%
PD FS Cor	2:41	1:06	1:35	59%
T2 FS Cor	2:41	1:46	:55	34%
PD FS Sag	3:32	1:06	2:26	69%
T1 Sag	3:23	1:03	2:20	69%
T1 FS Sag	3:14	1:46	1:28	45%
Lumbar	18:28	10:58	7:30	41%
T2 Sag	3:28	1:39	1:49	52%
T1 Sag	2:34	0:49	1:45	68%
T2 IR Sag	4:07	1:54	2:13	54%
T2 Ax	4:19	3:20	0:59	23%
T1 Ax	4:00	3:16	0:44	18%

Table 1 Data was acquired from the MAGNETOM Altea 1.5T System. The standard protocol was used as the "Pre Deep Resolve" exam and the "Post Deep Resolve" data incorporated the combination of Deep Resolve Gain & Sharp at similar protocol parameters. The decreased exam times as well as exam time percentage of reduction are included.

	1.5T Altea System	1.5T Altea with Deep Resolve	Estimated net increase with Deep Resolve Upgrade ¹
Average patients/day (7a - 7p)	12	17	5 more patients per day
Example reimbursement/patient²	\$450	\$450	\$450
Revenue/day (7a - 7p)	\$5,400	\$7,650	\$2,250
Revenue/month (24 days)	\$129,600	\$183,600	\$54,000
Revenue/year (288 days)	\$1,555,200	\$2,203,200	\$648,000

Table 2 Outlines potential increase in revenue with Siemens Deep Resolve if the decreased exam time allowed a customer to add 5 more patients per day.

[Learn more](#)



Siemens Healthineers Deep Resolve Boost

Radically shortens scan times without compromising image resolution

[siemens-healthineers.us](https://www.siemens-healthineers.us)

Deep Resolve Boost Competitive Scan Time Comparison

University Site Exam Type	Competitor		Siemens Healthineers		
	Standard Exam (min)	AI Based Accelerated Exam (min)	Deep Resolve Boost (min)	Improvement vs. Standard Competitor Exam	Improvement vs. AI based Accelerated Exam
Ankle	11:00	7:00	4:55	55%	30%
Knee	11:00	8:00	5:03	54%	37%
Shoulder	13:00	9:00	6:50	47%	24%
Hand/Wrist	19:00	13:00	7:25	61%	43%
C-Spine	15:00	11:00	4:44	68%	57%
T-Spine	12:00	8:00	3:59	67%	50%
L-Spine	15:00	9:00	4:28	70%	50%
Prostate	30:00	23:00	10:17	66%	55%
F. Pelvis	25:00	17:00	9:12	63%	46%



[Learn more](#)





Siemens Healthineers AG (listed in Frankfurt, Germany: SHL) pioneers breakthroughs in healthcare. For everyone. Everywhere. As a leading medical technology company headquartered in Erlangen, Germany, Siemens Healthineers and its regional companies is continuously developing its product and service portfolio, with AI-supported applications and digital offerings that play an increasingly important role in the next generation of medical technology. These new applications will enhance the company's foundation in in-vitro diagnostics, image-guided therapy, in-vivo diagnostics, and innovative cancer care.

Siemens Healthineers also provides a range of services and solutions to enhance healthcare providers' ability to provide high-quality, efficient care. In fiscal 2021, which ended on September 30, 2021, Siemens Healthineers, which has approximately 66,000 employees worldwide, generated revenue of €18.0 billion and adjusted EBIT of €3.1 billion.

Further information is available at www.siemens-healthineers.com.

The outcomes and statements provided by customers of Siemens Healthineers are unique to each customer's setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, and level of service/technology adoption), there can be no guarantee that others will achieve the same results.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features, which do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. For the most current information, please contact your local sales representative from Siemens Healthineers.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
siemens-healthineers.com

USA

Siemens Medical Solutions USA, Inc.
Healthcare
40 Liberty Boulevard
Malvern, PA 19355-9998, USA
siemens-healthineers.us